

# ACUVUE® Brand Contact Lenses for ASTIGMATISM Lens Calculator

- ✓ **9 out of 10** ACUVUE® Brand Contact Lenses for ASTIGMATISM lenses are **successfully fitted at first attempt**, with **minimal rotation from zero position**<sup>1-3</sup>
- ✓ **9 out of 10** Eye Care Professionals rated ACUVUE® Brand Contact Lenses for ASTIGMATISM **easy to fit**<sup>1</sup>

| Spectacle Spherical Power (DS) | ACUVUE® for ASTIGMATISM (sphere component in DS) | Spectacle Cylinder Power (DC) |           |           |           |           |           |           |           |           |       |       |       |       |       |
|--------------------------------|--|-------------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-------|-------|-------|-------|-------|
|                                |  | -0.75                         | -1.00     | -1.25     | -1.50     | -1.75     | -2.00     | -2.25     | -2.50     | -2.75     | -3.00 | -3.25 | -3.50 | -3.75 | -4.00 |
| +5.50                          | +6.00  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +5.25                          | +5.50  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +5.00                          | +5.25  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +4.75                          | +5.00  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +4.50                          | +4.75  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +4.25                          | +4.50  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +4.00                          | +4.25  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +3.75                          | +4.00  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +3.50                          | +3.50  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +3.25                          | +3.25  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +3.00                          | +3.00  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +2.75                          | +2.75  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +2.50                          | +2.50  |                               |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +2.25                          | +2.25  | cyl -0.75                     |           |           |           |           |           |           |           |           |       |       |       |       |       |
| +2.00                          | +2.00  |                               | cyl -1.25 |           |           |           |           |           |           |           |       |       |       |       |       |
| +1.75                          | +1.75  |                               |           | cyl -1.75 |           |           |           |           |           |           |       |       |       |       |       |
| +1.50                          | +1.50  |                               |           |           | cyl -2.25 |           |           |           |           |           |       |       |       |       |       |
| +1.25                          | +1.25  |                               |           |           |           | cyl -2.75 |           |           |           |           |       |       |       |       |       |
| +1.00                          | +1.00  |                               |           |           |           |           | cyl -2.75 |           |           |           |       |       |       |       |       |
| +0.75                          | +0.75  |                               |           |           |           |           |           | cyl -2.75 |           |           |       |       |       |       |       |
| +0.50                          | +0.50  |                               |           |           |           |           |           |           | cyl -2.75 |           |       |       |       |       |       |
| +0.25                          | +0.25  |                               |           |           |           |           |           |           |           | cyl -2.75 |       |       |       |       |       |

## Instructions

To select the appropriate ACUVUE® Brand Contact Lenses for ASTIGMATISM trial lens, cross-reference the patient's spectacle spherical power (in the far left-hand column with their spectacle cylinder power (along the top). For oblique cylinder axes, please refer to parameter chart on next page to check on availability.



| Spectacle Spherical Power (DS) | ACUVUE® for ASTIGMATISM (sphere component in DS) | Spectacle Cylinder Power (DC) |           |           |           |           |           |           |           |           |           |           |           |           |           |
|--------------------------------|--|-------------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|                                |  | -0.75                         | -1.00     | -1.25     | -1.50     | -1.75     | -2.00     | -2.25     | -2.50     | -2.75     | -3.00     | -3.25     | -3.50     | -3.75     | -4.00     |
| Plano                          | Plano  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -0.25                          | -0.25  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -0.50                          | -0.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -0.75                          | -0.75  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -1.00                          | -1.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -1.25                          | -1.25  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -1.50                          | -1.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -1.75                          | -1.75  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -2.00                          | -2.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -2.25                          | -2.25  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -2.50                          | -2.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -2.75                          | -2.75  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -3.00                          | -3.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -3.25                          | -3.25  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -3.50                          | -3.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -3.75                          | -3.75  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -4.00                          | -4.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -4.25                          | -4.00  | cyl -0.75                     |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -4.50                          | -4.25  |                               | cyl -1.25 |           |           |           |           |           |           |           |           |           |           |           |           |
| -4.75                          | -4.50  |                               |           | cyl -1.75 |           |           |           |           |           |           |           |           |           |           |           |
| -5.00                          | -4.75  |                               |           |           | cyl -2.25 |           |           |           |           |           |           |           |           |           |           |
| -5.25                          | -5.00  |                               |           |           |           | cyl -2.75 |           |           |           |           |           |           |           |           |           |
| -5.50                          | -5.00  |                               |           |           |           |           | cyl -2.75 |           |           |           |           |           |           |           |           |
| -5.75                          | -5.25  |                               |           |           |           |           |           | cyl -2.75 |           |           |           |           |           |           |           |
| -6.00                          | -5.50  |                               |           |           |           |           |           |           | cyl -2.75 |           |           |           |           |           |           |
| -6.25                          | -5.75  |                               |           |           |           |           |           |           |           | cyl -2.75 |           |           |           |           |           |
| -6.50                          | -6.00  |                               |           |           |           |           |           |           |           |           | cyl -2.75 |           |           |           |           |
| -6.75                          | -6.00  |                               |           |           |           |           |           |           |           |           |           | cyl -2.75 |           |           |           |
| -7.00                          | -6.50  |                               |           |           |           |           |           |           |           |           |           |           | cyl -2.75 |           |           |
| -7.25                          | -6.50  |                               |           |           |           |           |           |           |           |           |           |           |           | cyl -2.75 |           |
| -7.50                          | -6.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           | cyl -2.75 |
| -7.75                          | -7.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -8.00                          | -7.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -8.25                          | -7.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -8.50                          | -7.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -8.75                          | -8.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -9.00                          | -8.00  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -9.50                          | -8.50  |                               |           |           |           |           |           |           |           |           |           |           |           |           |           |
| -10.00                         | -9.00  | Try a spherical lens          |           |           |           |           |           |           |           |           |           |           |           |           |           |

Contact lens powers are calculated for back vertex distance (BVD) = 12.0 mm

(1) ACUVUE OASYS® for ASTIGMATISM (AOFA) (2) 1-DAY ACUVUE MOIST for ASTIGMATISM (1DAMFA) (3) ACUVUE OASYS® for ASTIGMATISM (AOFA) (4) ACUVUE VITA for ASTIGMATISM (AVFA)

