

Built on the ACUVUE® MOIST Platform—
the #1 prescribed daily disposable brand around the globe

Fit Success Tips


1-DAY ACUVUE® MOIST Brand MULTIFOCAL Contact Lenses:
Allowing You to Continue Meeting Your Patients' Eye Care Needs as They Age

1-DAY ACUVUE®
MOIST
BRAND CONTACT LENSES
MULTIFOCAL

ACUVUE® Brand has advanced the technology of multifocal lenses with 1-DAY ACUVUE® MOIST Brand MULTIFOCAL

EYE-INSPIRED™ Design:

a unique aspheric centre-near design for the presbyopic eye



A consistent and predictable fit helps you continue to address your patients' needs as they age
















Finding the best lens for the presbyopic eye

94% of patients were **successfully fit** with 2 pairs of lenses or less when using the Fitting Tips.*

Patient preparation:

1. Always start with a new subjective refraction, push the maximum plus for distance, and determine the spherical equivalent. Apply vertex distance correction if greater than +/-4.00D. Ensure that astigmatism is less than 1.00D.
2. Determine the dominant eye.
 - +1.00D blur test recommended
3. Determine the ADD based upon the patient's needs.
 - Focus on patient reading needs (desktop, tablet, mobile, etc)
4. Refer to the fit selection table based on the ADD to determine initial trial lenses.
 - Allow for 10 minutes of real-world exposure (outside of the exam room) before assessing visual performance
5. If your patient still has a specific near or distance need, refer to the enhancement tables.

To make your first lens the right lens, follow the selection table below:

For the best balance of distance and near vision			For patients requiring an enhancement for distance vision or near vision		
ADD	EYE	LENS SELECTION	ENHANCED DISTANCE VISION		ENHANCED NEAR VISION
+0.75D to +1.25D	Dominant Eye	 LOW	1-DAY ACUVUE® MOIST BRAND SPHERE	 LOW	 LOW
	Non-dominant Eye	 LOW			 LOW +
+1.50D to +1.75D	Dominant Eye	 MID	 LOW	 MID	
	Non-dominant Eye	 MID	 MID		 MID +
+2.00D to +2.50D	Dominant Eye	 MID	 MID	 MID	
	Non-dominant Eye	 HIGH	 MID +		 HIGH +

*: Add +0.25D to the distance power.



Allow for 10 minutes of real-world exposure (outside of the exam room) before assessing visual performance.

*Up to 4 lenses total.



1-DAY ACUVUE® MOIST Brand MULTIFOCAL

Product Information

LENS DETAILS

Material	etafilcon A
Diameter	14.3 mm
Base curve	8.4 mm
Wetting agent	Embedded PVP/LACREON® Technology
UV blocker*†‡	Class 2§
Water content	58%
Dk¶	28.0
Centre thickness	0.084 mm for -3.00D (varies with power)
Visibility tint	Yes
Inside-out mark	Yes
Powers	-9.00D to +6.00D, 0.25D steps
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D

PVP=polyvinylpyrrolidone.

† Helps protect against transmission of harmful UV radiation to the cornea and into the eye.

* **WARNING:** UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. You should continue to use UV-absorbing eyewear as directed. **NOTE:** Long-term exposure to UV radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of outdoor activities). UV-blocking contact lenses help provide protection against harmful UV radiation. However, clinical studies have not been done to demonstrate that wearing UV-blocking contact lenses reduces the risk of developing cataracts or other eye disorders.

‡ The UV blocking for these lenses averages 97% in the UVB range of 280 nm to 315 nm and 82% in the UVA range of 316 nm to 380 nm for the entire power range.

§ Class 2: Must transmit less than 50% (<50%) of UVA and 5% (<5%) of UVB.

¶ Units: $\times 10^{-11}$ (cm²/sec) (ml O₂/ml x mm Hg) @35°C Fatt (boundary corrected, non-edge corrected).

ACUVUE® Brand Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care, Division of Johnson & Johnson, Inc., by calling 1-800-267-5098 or by visiting www.acuvueprofessional.ca

ACUVUE®, 1-DAY ACUVUE® MOIST, EYE-INSPIRED™, INTUISIGHT™ and LACREON® are trademarks of Johnson & Johnson, Inc.
©Johnson & Johnson, Inc. 2015 April 2015 M-04-15-06-CE

