

# Now with **48%** more parameters to help cover even more of your astigmatic patients!<sup>1</sup>

**1-DAY ACUVUE®**  
**MOIST**  
 BRAND CONTACT LENSES  
**FOR ASTIGMATISM**

		LOW MINUS POWERS				HIGH MINUS POWERS			PLUS POWERS		
		Plano to -6.00 in 0.25D steps				-6.50 to -9.00 in 0.50D steps			+0.25 to +4.00 in 0.25D steps		
		-0.75	-1.25	-1.75	-2.25	-0.75	-1.25	-1.75	-0.75	-1.25	-1.75
OBLIQUES	30°	●	●	●							
	40°	●	●	●							
	50°	●	●	●							
	60°	●	●	●							
AGAINST THE RULE	70°	●	●	●	●	●	●	●	●	●	●
	80°	●	●	●	●	●	●	●	●	●	●
	90°	●	●	●	●	●	●	●	●	●	●
	100°	●	●	●	●	●	●	●	●	●	●
	110°	●	●	●	●	●	●	●	●	●	●
OBLIQUES	120°	●	●	●							
	130°	●	●	●							
	140°	●	●	●							
	150°	●	●	●							
WITH THE RULE	160°	●	●	●	●	●	●	●	●	●	●
	170°	●	●	●	●	●	●	●	●	●	●
	180°	●	●	●	●	●	●	●	●	●	●
	10°	●	●	●	●	●	●	●	●	●	●
	20°	●	●	●	●	●	●	●	●	●	●

● Current parameters    ● Expanded parameters

**2260**  
PARAMETERS

Covering  
**87%**  
 of astigmatic eyes<sup>1</sup>



**ACUVUE® Brand 1-Day Contact Lenses for ASTIGMATISM have the**

**LARGEST RANGE OF PARAMETERS**

**for astigmats among daily disposable contact lenses<sup>1</sup>**

1. JVC data on file July 2017. Proportion astigmats accommodated with ACUVUE® Brand Contact Lenses for ASTIGMATISM, and ACUVUE OASYS®, ACUVUE OASYS® 1-Day, ACUVUE® VITA & 1-DAY ACUVUE® MOIST Brand Contact Lenses (sphere & toric), and prevalence potential soft CL wearers requiring astigmatic correction.

ACUVUE® Brand Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care division of Johnson & Johnson Inc., by calling 1-800-267-5098 or by visiting acuvueprofessional.ca.

ACUVUE®, 1-DAY ACUVUE® MOIST and LACREON® are trademarks of Johnson & Johnson Inc. © Johnson & Johnson Inc. 2017 October 2017 MA-11-17-02-CE