

# ACUVUE® OASYS 1-Day with HydraLuxe® Technology

Unsurpassed coverage<sup>†</sup> of stock silicone hydrogel daily sphere and toric power combinations<sup>1</sup>

		LOW MINUS 0.00D to -6.00D (in 0.25D steps)				HIGH MINUS -6.50D to -9.00D (in 0.50D steps)			PLUS +0.25D to +4.00D (in 0.25D steps)		
		-0.75	-1.25	-1.75	-2.25	-0.75	-1.25	-1.75	-0.75	-1.25	-1.75
OBLIQUES	30°	✓	✓	✓							
	40°	✓	✓	✓							
	50°	✓	✓	✓							
	60°	✓	✓	✓							
AGAINST THE RULE	70°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	80°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	90°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	100°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	110°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
OBLIQUES	120°	✓	✓	✓							
	130°	✓	✓	✓							
	140°	✓	✓	✓							
	150°	✓	✓	✓							
WITH THE RULE	160°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	170°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	180°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	10°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓



**2260**  
TORIC  
PARAMETERS

## Extensive range of parameters

Unsurpassed coverage<sup>†</sup> available from stock lenses in Daily Disposable contact lenses to meet the needs of your patients.<sup>1</sup>



Around the clock coverage where you need it most— in low minus, plano to -6.00D, through -1.75DC\*.

NEARLY  
**96%**  
COVERAGE FOR  
SPHERICAL &  
ASTIGMATIC EYES

ACUVUE® OASYS Contact Lenses 1-Day with HydraLuxe® Technology family provide coverage for nearly 96% spherical and astigmatic eyes.<sup>1</sup>

<sup>†</sup> Unsurpassed coverage claim was demonstrated using a non-inferiority test with a difference margin of 5% and 95% confidence level.  
\* JUV data on file. Astigmatic coverage in Low Minus (-0.75, -1.25, -1.75DC across all axes in 10° increments) accommodated with ACUVUE® OASYS 1-Day with HydraLuxe® Technology for ASTIGMATISM.

<sup>1</sup> JUV Data on File 2020. Proportion Astigmats Accommodated with ACUVUE® Brand Contact Lenses for ASTIGMATISM.

**Important Safety Information:** ACUVUE® Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care division of Johnson & Johnson Inc., by calling 1-800-267-5098 or by visiting [jvisionpro.com](http://jvisionpro.com).

© Johnson & Johnson Vision Care division of Johnson & Johnson Inc. 2020 PP2020MLT6654 June 2020

