

# ACUVUE® OASYS

The widest parameter coverage for spherical and astigmatic eyes<sup>1</sup>

		LOW MINUS					HIGH MINUS					LOW PLUS					HIGH PLUS				
		Plano to -6.00D (in 0.25D steps)					-6.50D to -9.00D (in 0.50D steps)					+0.25D to +4.00D (in 0.25D steps)					+4.25D to +6.00D (in 0.25D steps)				
		-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75
OBLIQUES	30°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	40°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	50°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	60°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
AGAINST THE RULE	70°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	80°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	90°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	100°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	110°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
OBLIQUES	120°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	130°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	140°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	150°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
WITH THE RULE	160°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	170°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	180°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	10°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓



**4950**  
TORIC  
PARAMETERS

Accommodating nearly 99% of your spherical and astigmatic patients.<sup>1</sup>



Around the clock coverage where you need it most — in low minus, high minus and plus powers\* including 2.75 cyl

Confidently fit your patients with up to -3.50DC cylinder for ACUVUE® OASYS with ASTIGMATISM<sup>1,2</sup>

NEARLY  
**99%**  
COVERAGE FOR  
SPHERICAL &  
ASTIGMATIC EYES

ACUVUE® OASYS Contact Lenses family has more in stock toric & spheric parameters than any other brand.<sup>1</sup>



\* Available without a special order.

<sup>1</sup> It is the responsibility of the ECP to make the final decision on fit and acceptability, as there may be residual astigmatism.

<sup>2</sup> J.J.V. Data on File 2020. Proportion Astigmats Accommodated with ACUVUE® Brand Contact Lenses for ASTIGMATISM.

<sup>3</sup> J.J.V. Data on File 2020. ACUVUE® OASYS for ASTIGMATISM for Refractive Astigmatism of up to 3.50DC.

**Important Safety Information:** ACUVUE® Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop.

Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care division of Johnson & Johnson Inc., by calling 1-800-267-5098 or by visiting [invisionpro.ca](http://invisionpro.ca).

© Johnson & Johnson Vision Care division of Johnson & Johnson Inc. 2020 PP2020MLT6654 June 2020