

Now, the ACUVUE® VITA Brand Contact Lens family covers nearly 96% of spherical and astigmatic eyes¹

		Low Minus (25) 0.00D to -6.00D (in 0.25D steps)				High Minus (6) -6.50D to -9.00D (in 0.50D steps)			Plus (16) +0.25D to +4.00D (in 0.25D steps)		
		-0.75	-1.25	-1.75	-2.25	-0.75	-1.25	-1.75	-0.75	-1.25	-1.75
OBLIQUES	30°	✓	✓	✓							
	40°	✓	✓	✓							
	50°	✓	✓	✓							
	60°	✓	✓	✓							
AGAINST THE RULE	70°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	80°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	90°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	100°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	110°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
OBLIQUES	120°	✓	✓	✓							
	130°	✓	✓	✓							
	140°	✓	✓	✓							
	150°	✓	✓	✓							
WITH THE RULE	160°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	170°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	180°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	10°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20°	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓



2260
PARAMETERS

87% coverage for patients with astigmatism¹

NEARLY

96%

COVERAGE FOR SPHERICAL & ASTIGMATIC EYES¹

ACUVUE® VITA Brand Contact Lenses with HydraMax™ Technology family provide coverage for nearly 96% spherical and astigmatic eyes.¹

NEW ACUVUE®
vita
BRAND CONTACT LENSES
ASTIGMATISM

1. JJVC data on file 2017. Proportion astigmats accommodated with ACUVUE® Brand Contact Lenses for ASTIGMATISM, and ACUVUE OASYS®, ACUVUE OASYS® 1-Day, ACUVUE® VITA & 1-DAY ACUVUE® MOIST Brand Contact Lenses (sphere & toric), and prevalence potential soft CL wearers requiring astigmatic correction.

ACUVUE® VITA Brand Contact Lenses are indicated for vision correction as a daily wear lens with one-month recommended replacement. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care division of Johnson & Johnson Inc., by calling 1-800-267-5098 or by visiting acuvueprofessional.ca.

ACUVUE®, ACUVUE® VITA and HydraMax™ are trademarks of Johnson & Johnson Inc. © Johnson & Johnson Inc. 2017 August 2017 AV-08-17-04-CE